EzRay Air™ Portable

vatech EzRay Air™ Portable User manual

EzRay Air™ Portable

User Manual

Model : VEX-P250
Version : 1.45

- ENGLISH
- ESPAÑOL
- 日本語
- 한국어

For the Worldwide region

innovation **i**nside

"i" stands for 'innovation', one of the core values of VATECH, which aims to expand accessibility of medical solutions to more people.

ENGLISH

Notice.....	i
Conventions Used in this Manual	ii
1. System Overview	1
1.1 Indications for Use	1
1.2 Principles of Operation.....	1
1.3 Components.....	2
1.4 General View of the Equipment	3
2. Safety Instructions.....	9
2.1 General Safety Guidelines	9
2.2 Warnings and Safety Instructions	10
3. Operation.....	13
3.1 Power On/Off	13
3.2 Operation Mode	14
3.3 Positioning.....	19
3.4 Exposure	28
3.4.1 Remote Exposure Switch.....	29
3.5 Using the Battery	30
4. Cleaning and Maintenance	37
4.1 Cleaning	37
4.2 Maintenance.....	38
4.2.1 Maintenance Task Checklist.....	39
Appendix.....	41
A.1 Combining the Holder and the devices using tripod	41
A.2 Using the rotating rectangular cover	44

This page intentionally left blank.

Notice

The **VEX-P250** is a portable X-ray system.

This manual contains descriptions, operational instructions, imaging procedures for the **VEX-P250** dental X-ray system. It is recommended that you thoroughly familiarize yourself with this manual to make the most effective use of this equipment. Observe all cautions, safety messages and warnings that appear in this manual.

Due to continuous technological improvements, the manual may not contain the most updated information.

The illustrations/photos of the equipment in this manual are only for illustration purposes. The actual equipment may differ. For further information not covered in this manual, please contact us at:

VATECH Co., Ltd.

Phone: (+82) 1588 9510

E-mail: gcs@vatech.co.kr

This document is originally written in English.

The **VEX-P250** is referred to as Equipment in this manual.

Manual Name: EzRay Air Portable (Model: VEX-P250) User Manual

Version: 1.45

Publication Date: 2020-01

Conventions Used in this Manual

The following symbols are used throughout this manual. Make sure that you fully understand each symbol and follow the instructions which accompany it.

To prevent personal injury and damage to the equipment, please observe all warnings and safety information included in this document.

	<p>Indicates that a specific hazard is known to exist which through inappropriate conditions or actions may cause:</p> <ul style="list-style-type: none"> • Severe personal injury (to the operator and patient) • Substantial property damage.
	<p>Indicates that a potential hazard may exist which through inappropriate conditions or actions will or can cause:</p> <ul style="list-style-type: none"> • Minor injury • Property damage.
	<p>Indicates that a potential problem may exist which through inappropriate conditions or actions can cause:</p> <ul style="list-style-type: none"> • Property damage.
	<p>Indicates precautions or recommendations that should be used in the operation of the system, specifically:</p> <ul style="list-style-type: none"> • Using this Manual • Notes to emphasize or clarify a point.

1. System Overview

The VEX-P250, a portable dental X-ray system, operates on 21.6 V_{DC} supplied by a rechargeable Li-ion polymer battery pack. The portable X-ray system is an X-ray generating device that is mainly designed for dental examination (teeth and jaw). The portable X-ray system is composed of an X-ray generating part with an X-ray tube including a device controller, a power controller, a user interface, a beam limiting part, a backscatter shield, and an optional remote exposure switch. The VEX-P250 is designed to diagnose teeth and jaw through X-ray exposure using intra-oral image receptors.

1.1 Indications for Use


The device is a portable dental X-ray system intended for use by a radiographer, diagnostic radiographers, medical radiation technologists or healthcare professionals who specializes in the imaging of human anatomy for the diagnosis and treatment of pathology, for both adult and pediatric subjects for producing diagnostic dental radiographs for treatment of diseases of the teeth, jaw, and other oral structures using intra-oral image receptors.

1.2 Principles of Operation

X-rays are emitted when a high voltage is supplied to the X-ray tube assembly which frees electrons from the cathode. They hit anode to produce X-rays. The equipment acquires images by emitting X-rays continuously on the human tooth.


1.3 Components

No.	Item	Qty.
1	Main Body	1
2	Hand/Neck Strap	1
3	Power cord	1
4	Battery Charger	1


1.4 General View of the Equipment


Main Body


No.	Item	Description
1	Backscatter Shield	Shields from the backscattering radiation.
2	X-ray Beam Limiting Device	Limits the X-ray exposure area. Default type: Round Cone + Round Cover (FOV: Ø 6 cm)
3	X-ray Exposure Button	Press the button for X-ray exposure.
4	Handle	Grip the handle securely when using the system.
5	Power Button	Power On/Off button
6	Battery	Rechargeable Lithium-ion battery
7	Adaptor Connector	Connect the charging adaptor.
8	Remote X-ray Exposure Switch Port	Connect the X-ray exposure cable switch. Alternatively, can be used as a service port
9	Strap Loop	Connect the strap.
10	Control Panel	Display for the X-ray exposure settings and operation conditions
11	X-ray Generator	Includes the X-ray tube and the high-voltage generator.

Control Panel

No.	Item	Description
1		Tube Voltage/Current Indicator Indicates the tube voltage and tube current of the system.
2		Angle/Time Display Displays the X-ray exposure time, error code, cooling time and exposure angle.
3		Adult/Child Selection Indicates a patient type (adult or child).
4		Remaining Battery Indicator Indicates remaining battery level.
		Battery Charging Indicator 1) Indicates the battery charger is connected to the device. 2) Indicates the battery needs to be charged when it flickers .
5		SET Button Resets the X-ray exposure angle.
6		Speaker Sound alarm for the X-ray exposure

No.	Item	Description
7		<p>Tooth Type Selection</p> <p>Selects the tooth type.</p>
8		<p>X-ray Exposure Indicator</p> <p>Indicates the X-ray exposure status. (Green: Ready / Yellow: X-ray On)</p>
9		<p>Jog Dial</p> <p>Turn the jog dial left (-) or right (+) to select an X-ray exposure setting, press the jog dial to confirm the operating setting.</p>
10		<p>AUTO Button</p> <p>Selects a tooth and exposure time automatically based on the X-ray exposure angle.</p>

Available Option Items

No.	Figure	Option name
1		Backscatter Shield
2		Rectangular Cover 2x3 (3x2) FOV: 2x3 cm, 3x2 cm (This cover can be used as both 2x3 and 3x2.)
3		Rectangular Cover 4x3 (3x4) FOV: 4x3 cm, 3x4 cm (This cover can be used as both 4x3 and 3x4.)
4		Remote Exposure Switch
5		Rotating rectangular cover 4x3 (3x4) FOV: 4x3cm, 3x4cm (This adaptor can be used as both 4x3 and 3x4.)
6		Rotating rectangular cover 2x3 (3x2) FOV: 2x3cm, 3x2cm (This adaptor can be used as both 2x3 and 3x2.)
7		Base Holder*

No.	Figure	Option name
<p>* In case a tripod is used with the Base Holder, refer to the specifications below.</p> <ul style="list-style-type: none">• Fixing bolt size: 3/8 inch• Maximum supportable weight: about 5 kg• Minimum height: > 130 cm• More than 3 columns are required.• When using a tripod with 3 columns, make sure to have space of at least 1 m in width on the bottom.		

NOTICE

Refer to Appendix A.2 'Using the rotating rectangular cover' for instructions for use.

This page intentionally left blank.

2. Safety Instructions

2.1 General Safety Guidelines

- Mode of operation: Continuous operation with cyclic loading 60 times the exposure time wait before the next exposure can begin.
- This equipment is designed and manufactured to ensure maximum safety of operation. Operate and maintain it in strict compliance with the safety precautions and operating instructions contained in this manual.
- Only a legally qualified person, i.e., a radiographer, diagnostic radiographers, medical radiation technologists or healthcare professionals who specialize in the imaging of human anatomy for the diagnosis and treatment of pathology is allowed to operate this equipment in a controlled environment.
- Observe all local fire regulations. Always keep a fire extinguisher near the equipment.
- The equipment must be installed, maintained, and serviced by qualified service personnel according to the procedures and preventive maintenance schedules. Users can perform an only battery replacement.
- Ensure that the on/off switch is set to off when the equipment is not in use.
- Always disconnect the power supply before cleaning the equipment.
- DO NOT keep the equipment or its parts in a humid place or a liquid substance.
- Avoid placing the equipment near chemical storage and gas-filled storage facilities.
- The backscatter shield is provided as an optional item. If it is damaged or defective, contact your Service Representative for replacement.
- The backscatter shield (option) protects users from backscatter radiation that they might be exposed during X-ray exposure. Operating the equipment with a backscatter shield can expose you to less radiation when operating without a backscatter shield.

2.2 Warnings and Safety Instructions

This X-ray unit may be dangerous to patient and operator unless safe exposure factors, operating instructions, and maintenance schedules are observed.

It is essential to read this user manual carefully and strictly abide by all warnings and cautions stated within it.

To avoid electric shock, this equipment must be connected to a mains power source with a protective earthing connection.

Since rules and regulations concerning radiation safety differ between countries, it is the responsibility of the owner and operator of this equipment to comply with all applicable rules and regulations concerning radiation safety and protection in their area.

- DO NOT open or remove the cover panels on this equipment.
- DO NOT expose this equipment to liquids, mists or sprays. Exposing this equipment to liquids may cause an electrical shock or otherwise damage the system.
- DO NOT use spray cleaners on this equipment, as this could cause a fire.
- DO NOT use this equipment in an environment that is susceptible to explosion.
- DO NOT place flammable materials near this equipment.
- DO NOT touch the patient while also touching the SIP/SOP connectors.
- DO NOT try to modify this equipment, including the wires or cables. Modifying this equipment may damage it beyond repair.
- it is recommended that the patient and the operator wear protective lead-lined aprons unless other Radiation Protection Protocols apply locally.
- Children and pregnant women must consult with a doctor before X-ray exposure.
- Grave dangers may occur from electromagnetic interference (i.e., noise) between other equipment in the area during specific examinations or medical treatment.
- Battery chargers must be located in an accessible area where they can be easily unplugged from the power source.

Battery Use

- Make sure to charge the battery in the external environment from the patient.
- Make sure to use the battery only provided or approved by VATECH. If non-standard or damaged batteries are used, there is a risk of fire and explosion.
- Make sure to use the battery charger only provided or approved by VATECH. Using an unauthorized charger may result in battery damage.
- DO NOT expose batteries to heat or fire. Avoid storage in direct sunlight.
- DO NOT short-circuit, crush, puncture, mutilate, or disassemble the battery.
- DO NOT store batteries haphazardly in a box or drawer where they may short-circuit each other or be short-circuited by other metal objects.
- Observe the plus (+) and minus (-) marks on the battery and equipment and ensure correct use.
- DO NOT remove a battery from its original packaging until required for use.
- DO NOT subject batteries to mechanical shock.
- In the event of a cell leaking, do not allow the liquid to come in contact with the skin or eyes. If the contact has been made, wash the affected area with copious amounts of water and seek medical advice.
- DO NOT make the battery wet or let it be in the water. Keep batteries clean and dry.
- Keep batteries away from children and pets.
- Seek medical advice immediately if a battery has been swallowed.
- DO NOT dispose of batteries with ordinary trash. Turn in discharged batteries to local supply or discard or recycle batteries according to your local government regulations.
- Batteries can be replaced by users.
- Make sure to turn off the device before replacing the battery.
- When charging the battery, the exposure function is locked.
- DO NOT leave a battery on a prolonged charge when not in use.
- If the equipment has not been used for long periods of time, it is recommended to charge the battery before use.
- After extended periods of storage, it may be necessary to charge and discharge the cells or batteries several times to obtain maximum performance.
- Be sure to turn off the equipment when not in use. This helps to ensure the life of the battery.

2. Safety Instructions

- Be sure to charge the battery frequently. This helps to ensure the life of the battery.
- If the equipment not in use has been turned on for long periods of time, the battery may be fully discharged.
 - Depending upon the battery discharge status, it takes about 1 day for charging the battery. If the device is not turned on after charging the battery for about 1 day, it indicates that the battery has been fully discharged. Contact your Service Representative for battery replacement.
 - DO NOT charge a fully discharged battery, as this may cause fire or explosion. Be sure to replace the battery (provided by VATECH).

Radiation Safety


When using the equipment, it is recommended that all users comply with the following radiation safety guidelines for the safety of the users and the patients.

- All users and patients should wear protective equipment, such as a lead apron, thyroid collar, et cetera
- This equipment should be operated in the area that is more than 6 feet away from other personnel, such as assistants or other patients. If they should stay closer than 6 feet, it is recommended that they wear a lead apron, thyroid collar, or stay behind a lead shield.
- Pregnant women should not be exposed to X-rays unless it is strictly necessary.
- All users should comply with the Radiation Protection Policies established by the government.
- When selecting a Position Indicating Device, it should be considered if the PID can be used with the backscatter shield attached at the outer end of the cone for maximum operator protection.

3. Operation

3.1 Power On/Off

1. Turn on the system referring to the following figure.


2. The following displays and indicators light up:
 - Current Angle/Time display
 - Tooth type selection display
 - Adult/Child selection display
 - Remaining battery indicator
 - X-ray exposure indicator
3. Make sure that at least one battery indicator light comes on.

Battery level 1

NOTICE

When the battery indicator has one flickering light, charge the battery immediately by using the battery charger.

For more information, refer to '3.5 Using the Battery'.

3.2 Operation Mode


This system can be operated with Manual Mode and Auto Mode, and you can set up the mode by using the **AUTO** button.

Manual Mode


1. To start the Manual Mode, check if the lamp under the **AUTO** button is turned off.

2. When the tooth type selection area flickers, turn the jog dial to select the tooth type. To see the Control Panel before and after selection, refer to the figures below.


Tooth Type

Symbol	Type
	Incisor
	Canine
	Molar/Premolar
	Bitewing

- After tooth type selection, a patient type should be selected. When the Adult/Child selection area flickers, turn the jog dial to select the patient type. To see the Control Panel after selection, refer to the figure below.

After patient type selection

Patient Type

Symbol	Type
	Adult
	Child

NOTICE

After the tooth type and patient type are selected, the exposure time is automatically displayed.

- If you want to change the exposure time, turn the jog dial to adjust the exposure time from 0.05 to 1.0 s. (increments: 0.01 s)


NOTICE

If you press the jog dial after adjusting the exposure time in Manual Mode, the exposure time is returned to the default setting.

To save the exposure time as default in Manual Mode, press and hold the jog dial for about 3 seconds.

Auto Mode

1. When the Auto Mode is turned on by pressing **AUTO**, the default angle is displayed as shown in the following figure.

NOTICE

To set the starting point (0°) during exposure, press **SET**.

2. Select the patient type by turning the jog dial.
3. Perform patient positioning and press **SET**.

To check how to perform patient positioning, refer to '3.3 Positioning' on *page 19*. The following figure is an example of the maxilla incisor.

- When tooth types are selected, the exposure angles are automatically displayed according to the tooth type. To check the default exposure angles, refer to the following table.

Exposure time is automatically set according to the patient and tooth type.

Tooth type	Angle of inclination
Incisor	Maxilla: +40° ~ +50°
	Mandible: -22° ~ -28°
Canine	Maxilla: +40° ~ +50°
	Mandible: -17° ~ -23°
Molar/Premolar	Maxilla: +25° ~ +35°
	Mandible: -2° ~ -8°
Bitewing	+3° ~ +12°

NOTICE

Since the angles of inclination for the maxilla incisor and canine are the same, the exposure time of the canine is applied to both the maxilla incisor and canine.

Refer to the following figure to see the angle for molar/premolar.


3.3 Positioning


Positioning the Patient

To obtain high-quality intra-oral radiography with maximum details, take extra care in all steps of the radiography process: positioning the patient and the X-ray imaging system; exposing the phosphor plate or the intra-oral sensor.

1. Place a protective lead apron on the patient's chest.
2. Have the patient sit on the chair with the sagittal plane vertical.
 - For radiography of the upper maxillary, the Frankfort plane must be horizontal.

- For radiography of the lower maxillary, the Occlusal plane must be horizontal.

3. Place the tube head cone on the area you want to take an image.

When you hold the device, it is recommended to grip the handle by one hand and place the other on the underside between the Exposure Button and the cone as shown in the following figure.


NOTICE

Depending on the imaging angles, exposure times vary. Since it is necessary to keep the patient with low X-ray doses and the user in the protected area, have the patient's head slightly tilted, and raise or lower the chin if needed. Please refer to '2.2 Warnings and Safety Instructions'.


For further information about the patient positioning and beam angle for each mode, refer to the following **Positioning Instructions**.

Positioning Instructions

Paralleling technique: The phosphor plate or sensor is placed in a holder that is used to align the phosphor plate or sensor parallel to the long axis of the teeth.

Bisected angle technique: The patient holds the phosphor plate or sensor in place with his/her finger. The X-ray beam is directed perpendicularly towards an imaginary line, which bisects the angle between the phosphor plate or sensor plane and the long axis of the tooth.

Position the tube head to the patient using the accepted standard positioning procedures.


Here are the specific angulations and directions for the tube head to take the best images of a particular tooth (i.e., **Bisected angle technique**).

Position the receptor carefully not to damage the soft tissue of the patient's intra-oral area.

- **Maxillary Incisor**


The x-ray beam is directed downward at 45°.

Teeth		Angle of inclination
Incisor	Maxilla	+45°

- **Mandibular Incisor**


The x-ray beam is directed upward at 25°.

Teeth		Angle of inclination
Incisor	Mandible	-25°

- **Maxillary Canine**


The x-ray beam is directed downward at 45°.

Teeth		Angle of inclination
Canine	Maxilla	+45°

- **Mandibular Canine**


The x-ray beam is directed upward at 20°.

Teeth		Angle of inclination
Canine	Mandible	-20°

- **Maxillary Molar and Premolar**


The x-ray beam is directed downward at 30°.

Teeth		Angle of inclination
Molar and Premolar	Maxilla	+30°

- **Mandibular Molar and Premolar**

The x-ray beam is directed upward at 5°.

Teeth		Angle of inclination
Molar and Premolar	Mandible	-5°

- **Bitewing**

For a bitewing exposure, the patient closes their teeth during exposure on the phosphor plate/sensor holder.

The x-ray beam is directed downward at $5^{\circ} \sim 8^{\circ}$.

Teeth	Angle of inclination
Bitewing exposure	$+5^{\circ} \sim +8^{\circ}$

Positioning the Imaging Sensor

Using the VEX-P250 dental X-ray system, you may create an X-ray image on different types of imaging receptors:

- Digital sensors
- Phosphor plate


The digital imaging sensor must be positioned appropriately to ensure image quality. (for information about the proper placement of the imaging sensor, please refer to 'Positioning Instructions'.)

- Failure to position the imaging sensor correctly can result in errors on the radiograph, such as distorted teeth and roots, elongation, magnification, and overlapping contacts.


NOTICE

The paralleling technique generally reduces the risk of such errors, but if you position the sensor improperly, angulation errors may occur (angulation of the sensor to the tooth itself).


- Failure to align the imaging sensor with the exit pattern of the X-ray beam can result in cone-cuts on the radiograph. The cone-cuts are bright areas that are shown on the radiograph when part of the radiograph is not exposed to radiation. Please refer to the following figure as an example of cone-cuts.


The following figure indicates how the cone-cut occurred by showing the position of the imaging sensor and the radiation area.

To ensure proper alignment between the imaging sensor and the X-ray beam, it is recommended to use a PID (Position Indicating Device) as shown in the following figure.

When using the PID, the exit pattern of the X-ray device should be aligned perpendicular to the target receptor as shown in the following figure.

**NOTICE**


Once the PID is appropriately aligned, instruct the patient not to move.

3.4 Exposure

IMPORTANT

The operator **MUST** instruct the patient to refrain from moving during the entire exposure.

1. Instruct the patient not to move.
2. Press the Exposure Button for exposure duration.

3. While X-ray is being exposed,
 - The X-ray Exposure Indicator lights up, and an audible sound is produced.
 - Keep pressing until the X-ray Exposure Indicator light goes out and the audible sound stops.

Green: Ready

Yellow: X-ray On

IMPORTANT


Hold the Exposure Button or switch as long as the acoustic signal can be heard. Otherwise, the exposure will be faulty, and there will be an error message on the Control Panel.

3.4.1 Remote Exposure Switch

The **Remote Exposure Switch** allows the operator to control image acquisition from outside of the X-ray room.

Press and hold the **Remote Exposure Switch** until the acquisition is completed. Premature release of the **Remote Exposure Switch** will abort image acquisition.

Pressing the **Remote Exposure Switch** activates the X-ray Exposure Indicator to turn yellow. This color indicates that the X-ray is being emitted.

IMPORTANT

The **Remote Exposure Switch** is detachable. Ensure that the **Remote Exposure Switch** cable is not detached out from the unit accidentally during the operation.

IMPORTANT

Keep vocal/visual contact with the patient during exposure. If any problem occurs during exposure, release the **Remote Exposure Switch** immediately.

3.5 Using the Battery

Battery level indicator with a residual quantity is shown on the left side of the Control Panel. When the battery indicator has one flickering light (level 1), charge the battery immediately. Refer to the battery levels as shown below.

Level 3

Level 2

Level 1

If the battery is not charged for about an hour, the battery voltage becomes low. At this level, all displays are turned off except for the error code **A.10** as shown in the following figure, so exposure cannot be performed. After the device is connected to the battery charger, and the battery level becomes 'Level 1', all functions are returned to normal operation.

IMPORTANT

System Status depending on Battery Levels

Item	System Status				
	Battery Level 3, 2	Battery Level 1		Low Battery	
When turning on the system	Operation	Normal	Normal	Normal	Not operated
	Battery Level Indicator	Normal	Normal	Flickers	Not displayed
	Battery Charging Indicator	Not displayed	Not displayed	Flickers	Not displayed
	Control Panel Brightness	Normal	Normal	Dark	Normal (Error code A.10 is only displayed)
When operating the system	Operation	Normal	Normal	Normal	Not operated
	Battery Level Indicator	Normal	Normal	Flickers	Not displayed
	Battery Charging Indicator	Not displayed	Not displayed	Flickers	Not displayed
	Control Panel Brightness	Normal	Normal	Normal	Normal (Error code A.10 is only displayed)

Make sure to charge the battery if the Battery Charging Indicator flickers (or the error code A.10 is displayed). If the device has been turned on for long periods of time with the error code A.10, the battery may be discharged.


DO NOT charge a fully discharged battery, as this may cause fire or explosion. Be sure to replace the battery (provided by VATECH).


When the device is connected to the battery charger, the Battery Charging Indicator is displayed at all times except that the battery is fully discharged.

Charging the Battery

1. Connect the battery charger to the battery charger connector as shown in the following figure.

1. When the battery charger is connected, the battery charging LED indicator light comes on. Charge the battery until all the three LED indicators are filled up.

NOTICE

Usually, it takes about 3 hours to charge the battery after a complete discharge fully.

2. When the battery charge is completed, remove the battery charger from the device.

NOTICE

You cannot perform an exposure while the battery charger is connected to the device.

CAUTION

Please use only the latest charger supplied with this device since charging specifications is differ from the old charger. Charging this device using an older charger may damage the device.

Battery Use Cycle

The battery is a consumable part. It is expected to degrade gradually, so it should be recharged more frequently. When the battery duration decreases to half or less than half compared to when the battery was new, contact your Service Representative to get a new battery.

To check how to replace the battery, see 'Battery Replacement' on page 33.


Battery Replacement

NOTICE


Users can replace batteries.

To replace the battery, contact your Service Representative to get a battery kit (including a new battery and a Phillips screwdriver).


1. Workforce
One person, 3 minutes
2. Tools required
Phillips screwdriver (size: T20)
3. Removal Procedure
 - 1) Using a Phillips screwdriver, unscrew the battery bay access door.

- 2) Lift the door and remove the battery.

- 3) Disconnect the battery cable from the device by pressing the battery cable connector as shown in the following figure.


DO NOT pull excessively on the battery cable.

- 4) Install the new battery in the reverse order of removal.

Sleep Mode

To reduce battery consumption, the Sleep Mode is started when the system is not in use for one minute.

1. When the Sleep Mode is started, the Control Panel becomes dark as shown in the following figure (right).

2. To return to normal operation, move the system or press any button (except for the X-ray Exposure Button).
3. If you do not use the system, turn it off to reduce battery consumption.

Power Off Mode

If the system (not in use) has been turned on for 30 minutes since the Sleep Mode was started, the Power Off Mode is started.

1. 'Off' is flickering on the Control Panel. If you need to be notified by an alarm when the Power Off Mode is started, change the setting.

2. To return to normal operation, move the system or press any button (except for the X-ray Exposure Button).
3. If you do not use the system, turn it off to reduce battery consumption.

Power Save Mode

If the system (not in use) has been turned on for 4 hours since the Power Off Mode was started, the Power Save Mode is started.

1. All displays are turned off on the Control Panel.
2. To return to normal operation in the Power Save Mode, you **MUST** turn off the system and turn it back on.

4. Cleaning and Maintenance

4.1 Cleaning

WARNING

Before cleaning the equipment, make sure to turn off the equipment.

- The equipment surfaces can be cleaned with a soft cloth dampened in an alcohol-based, non-corrosive cleaning solution. If necessary, wipe off surfaces with disinfectant.
- If necessary, wipe off surfaces with disinfectant.
- Please observe the hygiene instructions of the phosphor plate scanner manufacturer.

CAUTION

When cleaning the surfaces, make sure that the equipment is not connected to the battery charger.

CAUTION

- DO NOT expose the equipment to any liquids.
- DO NOT use spray cleaner or disinfectant directly into the equipment as this could cause a fire.

NOTICE

The soft cloth should be damp, but not dripping wet.

NOTICE

The cloths or wipes cannot be re-used.

4.2 Maintenance

VATECH requires periodic constancy tests to ensure image quality and safety for the patient and operator.

Only **VATECH** authorized technicians can perform inspection and service of this equipment. For technical assistance, contact the **VATECH** service center or your local **VATECH** representative.

Please observe the hygiene instructions of the phosphor plate scanner manufacturer.

Cautions and Notes

DO NOT keep the equipment or its parts in a humid place or near a liquid substance.

Avoid placing the equipment near chemical storage and gas-filled storage facilities.

When the equipment is not in use for a long time, fully charge the battery and remove it from the device before storage.

4.2.1 Maintenance Task Checklist

Always turn off the equipment before performing any maintenance.

Tasks	Period
Before the operation, ensure that the equipment is clean and ready for use.	Daily
After using the equipment, make sure that the equipment has been turned off.	Daily
Wipe the outer covers of the equipment with a dry cloth at the end of each day's operation.	Daily
 DO NOT use detergents or solvents to clean the outer covers of the equipment.	Daily
Ensure that the signal is audible and the X-ray emission light is visible when you make an exposure.	Daily
Ensure that the yellow (exposure) indicator light turns on when the Exposure Button is pressed.	Daily
Ensure that the battery charging LED indicator comes on when charging the battery.	Daily
Ensure that the battery level indicator displays at least two levels (Battery Level 2). For more information on the battery levels, see '3.5 Using the Battery' on <i>page 30</i> .	Daily
Ensure that all visible labels are intact and legible.	Monthly

If any defects are found, do not operate the equipment since it has to be handled by a qualified person. Contact your Service Representative.

This page intentionally left blank.

Appendix

A.1 Combining the Holder and the devices using a tripod

NOTICE

Portable x-rays are not available in some countries. To be used in these countries, Vatech provides the holder with the exposure switch as an option.

The cradle is packed/shipped in a separate box as shown in the image below.

Remove all the components in the box and follow the procedure below to attach the cradle. (Note: Tripod not included)

NOTICE

To remove the holder, reverse the mounting procedure.


1. Using a Phillips screwdriver, connect the EzRay Air Portable and Gender cable with the supplied screws.

2. Attach the holder to a tripod. Refer to the Tripod manual for how to attach to a tripod.

3. Attach the EzRay Air Portable unit to the holder with the tripod as shown in the image below.

4. Connect one end of the Gender cable to the EzRay Air Portable device.


5. Connect the other end of the Gender cable connected to the holder to the supplied exposure cable.


A.2 Using the rotating rectangular cover

Rotating rectangular cover rotates in 360 degrees.

Rotating rectangular cover consists of the Outer Ring and the Inner Ring.


The Outer Ring rotates in 90-degree increments.

IMPORTANT

The Inner Ring rotates in 10-degree increments.

1. Assemble the rotating rectangular cover to the Cone of the Main body.

NOTICE

For assembling the rotating rectangular cover to the Cone, turn the rotating rectangular cover 45 degrees and check the Click sound to make sure the assembly is complete.

2. Turn the Outer Ring of the rotating rectangular cover to adjust the angle roughly.
3. Turn the Inner Ring of the rotating rectangular cover to fine-tune the angle.

NOTICE

The disassembly method of the rotating rectangular cover is the same as the assembly method. Refer to the NOTE mentioned below stage 1.

Copyright by © 2016 VATECH Co., Ltd.

All rights reserved.

The documentation, brand name, and logo used in this manual are copyrighted.

No part of this manual may be reproduced, transmitted, or transcribed without the expressed written permission of the manufacturer.

We reserve the right to make any alterations that may be required due to technical improvement. For the most current information, contact your VATECH representative.

Manufactured by VATECH Co., Ltd.

Tel: (+82) 1588 9510

Email: gcs@vatech.co.kr

Website: www.vatech.co.kr

Head Quarters Address: 13, Samsung 1-ro 2-gil, Hwaseong-si, Gyeonggi-do, 18449, Korea

Factory Address: 13, Samsung 1-ro 2-gil, Hwaseong-si, Gyeonggi-do, 18449, Korea

Authorized EU Representative

Name: Vatech Dental Manufacturing Ltd.

Address: Suite 3, Ground Floor, Chancery House, St. Nicholas Way, Sutton, SM1 1JB UK

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Copyright © 2016 VATECH Co., Ltd.