

3D DENTAL —SYSTEMS—



3DCART

3D DENTAL —SYSTEMS—

3D Dental Systems d.o.o.

Delavska ulica 19 | 3000 Celje | Slovenija | EU

Phone: +386 40 634 802

Email: info@3dds.eu




www.3dds.eu

ABOUT US

Our young company is a team of highly driven people who develop intuitive solutions for dental practices and laboratories.

Our goal is to provide you with modern-looking and high-quality products for professional use.

All our products are designed and produced in Slovenia (EU). We use only the best possible material for every part, and all components we don't produce ourselves are heavily tested before we decide to use them in our products. With all of this in mind, you can be sure that all our products are of the highest quality.




THE DESIGN

We are proud to present to the **3D Cart**.

Modern clean lines combined with **quality materials** will enable this cart to **stand out** in your dental office.

3D Cart is designed to ensure that only power and intraoral scanner cables are visible, hiding others inside the cart. This allows your scanning setup to **look professional** and makes **every-day work easier**.



INTUITIVE

The scanner is positioned in a way that makes it **easily accessible** from your workspace, **yet safely hidden**, so you can't bump into it by mistake.



STRONG AND RELIABLE

The body is constructed from **6 mm thick aluminum**, making it strong and just enough heavy.

Casters have a 360° rotating design with lock function on the front two, allowing the 3D Cart to be mobile or fixed.

STORAGE FOR ACCESSORIES

The shelf inside the 3D Cart offers **additional storage** space for scanning tips and calibration tools.

There is sufficient space inside the cart for both PC and UPS, if needed.

FREEDOM OF CHOICE

You can use the 3D Cart **with a desktop or laptop computer** without the need to buy any additional add-ons. This gives you the flexibility in deciding which computer would be most suitable for your needs and gives you more options when you decide to upgrade your computer.

The 3D Cart can be used **with any intraoral scanner** that uses a PC or laptop (Medit, CareStream, 3Shape, Planmeca, Shining3D and more).

OPTIONAL ACCESSORIES

DISPLAY SUPPORT

Display support with basic monitor mount is always included with the 3D Cart. With additional mount you can add tilt option to your monitor.



MONITOR MOUNT

White monitor mount with tilt option.




POWER EXTENSION CORD


White spiral power extension cord.
Two available lengths:
- 2.5 m (3 sockets)
- 5.0 m (4 sockets)

SPECIFICATION AND DIMENSIONS

Base size	83 x 55 x 40 cm
Work surface	55 x 40 cm
PC cabinet size	50 x 29 cm
Weight	25 kg
Colour	White, powder coat
Display support	VESA monitor mount - up to 32"
Body	6mm aluminum
Doors	Plexiglass
Wheels	360° ; lock function on front two



FRONT VIEW



SIDE VIEW

3D DENTAL —SYSTEMS—

3D Dental Systems d.o.o.

Delavska ulica 19 | 3000 Celje | Slovenija | EU

Phone: +386 40 634 802

Email: info@3dds.eu

www.3dds.eu