



2019

codent
Let's smile together

APOLLO SERIES A6 AIR TURBINES

four sprays - dual coat titanium coating - ceramic bearings - 2-year warranty



MODEL A6L

Standard head - NSK® PTL
®QUICK coupling
Burs: 1.6 - Cellular glass optic
Rotation: 420,000 rpm
Head size: Ø12.4 x H12.9 mm



LED

MODEL A6N

Standard head - NSK®MACH
®QUICK coupling
Burs: 1.6 - No cellular glass optic
Rotation: 380,000 rpm
Head size: Ø12.4 x H12.9 mm



MODEL A6KL

Standard head - KaVo®
MULTIflex® coupling
Burs: 1.6 - Cellular glass optic
Rotation: 380,000 rpm
Head size: Ø12.4 x H12.9 mm



LED

MODEL A6K

Standard head - KaVo®
MULTIflex® coupling
Burs: 1.6 - No cellular glass optic
Rotation: 380,000 rpm
Head size: Ø12.4 x H12.9 mm



MODEL A6BL

Standard head - BIEN-AIR® UNIFIX® coupling
Burs: 1.6 - Cellular glass optic
Rotation: 380,000 rpm
Head size: Ø12.4 x H12.9 mm




LED



MODEL A6WL

Standard head - WH® ROTO QUICK® coupling
Burs: 1.6 - Cellular glass optic
Rotation: 380,000 rpm
Head size: Ø12.4 x H12.9 mm



LED



APOLLO SERIES A5 AIR TURBINES

four sprays - dual coat titanium coating - ceramic bearings - 2-year warranty



MODEL A5L

Small head - NSK® PTL
® QUICK coupling
Burs: 1.6 - Cellular glass optic
Rotation: 420,000 rpm
Head size: Ø11.2 x H12.5 mm

MODEL A5N

Small head - NSK® MACH
® QUICK coupling
Burs: 1.6 - No cellular glass optic
Rotation: 420,000 rpm
Head size: Ø11.2 x H12.5 mm

LED



MODEL A5KL

Small head - KaVo®
MULTIflex® coupling
Burs: 1.6 - Cellular glass optic
Rotation: 420,000 rpm
Head size: Ø11.2 x H12.5 mm

MODEL A5K

Small head - KaVo®
MULTIflex® coupling
Burs: 1.6 - No cellular glass optic
Rotation: 420,000 rpm
Head size: Ø11.2 x H12.5 mm


LED



MODEL A5BL

Small head - BIEN-AIR® UNIFIX® coupling
Burs: 1.6 - Cellular glass optic
Rotation: 420,000 rpm
Head size: Ø11.2 x H12.5 mm

LED



MODEL A5WL

Small head - WH® ROTO QUICK® coupling
Burs: 1.6 - Cellular glass optic
Rotation: 420,000 rpm
Head size: Ø11.2 x H12.5 mm

LED



APOLLO A6 - A5 SERIES AIR TURBINES 4 HOLES & B2

four sprays - dual coat titanium coating - ceramic bearings - 2-year warranty



MODEL A6M4

Standard head - 4- holes connection
Burs: 1.6 - No cellular glass optic
Rotation: 360 000 rpm ~ 420 000 rpm
Head size: \varnothing 12.5 x H12.6 mm



MODEL A6B2

Standard head - B2 - 2 holes connection
Burs: 1.6 - No cellular glass optic
Rotation: 360 000 rpm ~ 420 000 rpm
Head size: \varnothing 12.5 x H12.6 mm



MODEL A5M4

Small head - 4- holes connection
Burs: 1.6 - No cellular glass optic
Rotation: 380 000 rpm ~ 450 000 rpm
Head size: \varnothing 10.9 x H11.9 mm



MODEL A5B2

Small head - B2 - 2 holes connection
Burs: 1.6 - No cellular glass optic
Rotation: 380 000 rpm ~ 450 000 rpm
Head size: \varnothing 10.9 x H11.9 mm



APOLLO SERIES A45 AIR TURBINES

four sprays - dual coat titanium coating - ceramic bearings - 2-year warranty



LED

MODEL A45L

45° head - NSK® PTL® QUICK coupling
Burs: 1.6 - Cellular glass optic
Rotation: 450,000 rpm
Head size: Ø11 x H13.5 mm



LED

MODEL A45KL

45° head - KaVo® MULTIflex® coupling
Burs: 1.6 - Cellular glass optic
Rotation: 450,000 rpm
Head size: Ø11 x H13.5 mm

MODEL A45K

45° head - KaVo® MULTIflex® coupling
Burs: 1.6 - No cellular glass optic
Rotation: 450,000 rpm
Head size: Ø11 x H13.5 mm



MODEL A45 M4

45° head - 4-hole connection
Burs: 1.6 - No cellular glass optic
Rotation: 450,000 rpm
Head size : Ø11 x H13.5 mm





Q6CM- LED COUPLING

Compatible with NSK® PTL® instruments
6-holes - LED illumination

Q4M COUPLING

Compatible with NSK® MACH® instruments
4-holes - No light



Q6K- LED COUPLING

Compatible with KaVo® MULTiflex® instruments
6-holes - LED illumination



Q6KW- LED COUPLING


























Compatible with KaVo® MULTiflex® instruments
6-holes - LED illumination - Water controller



Q4K COUPLING

Compatible with KaVo® MULTiflex® instruments
4-holes - Without illumination

instruments compatible with couplings

	A6L	Q4M Q6CM-LED	
	A6N		
	A5L		
	A5N		
	A45L		
	A6KL	Q4K	
	A6K		
	A5KL	Q6K-LED	
	A5K		
	A45KL	Q6KW-LED	
	A45K		
	A6WL	W&H Coupling*	
	A6BL	Bien air Coupling*	
	APJ BA	Bien air Coupling*	
	APJ WH	W&H Coupling*	
	APJ KV	Q6K-LED Q6KW-LED Q4K	
	APJ NM	Q4M Q6CM-LED	
	APJ M4	4 HOLES	

* The Bien-Air and W&H couplings are cited as examples but not marketed by Codent-Europe GmbH



codent
Let's smile together